



Technical Education, Vocational and Entrepreneurship  
Training Authority (TEVETA)

**DIPLOMA IN SCIENCE LABORATORY TECHNOLOGY**

**YEAR I**

**Chemistry Techniques I**

**Record of Practical Assessment**

Learner`s Name:\_\_\_\_\_

Learner`s NRC no.:\_\_\_\_\_

Learner`s TEVETA No.:\_\_\_\_\_

Assessment Period:\_\_\_\_\_

Copyright

## **PREFACE**

The Technical Education, Vocational and Entrepreneurship Training Authority (TEVETA) is an institution created under the Technical Education, Vocational and Entrepreneurship Training Act Number 13 of 1998, as amended by the Technical Education, Vocational and Entrepreneurship Training (Amendment) Act Number 11 of 2005.

The Act among other things provides that TEVETA shall:

- (a) regulate and conduct national examinations and assessments relating to technical education, vocational and entrepreneurship training;
- (b) charge and collect fees in respect of examinations, assessments and other services provided by the Authority;
- (c) award certificates to persons who succeed in examinations and assessments undertaken under this Act
- (d) do all such things connected with or incidental to the functions of the Authority under this Act.

Through this mandate, the Assessment and Qualifications Division of TEVETA has developed Practical Assessment Tool Kits to enable learners achieve the competences that are congruent with the demand of the workplace tasks. These tool kits in part are also intended to ensure that similar conditions under which all students in TEVET are assessed and examined apply wherever the course is undertaken in Zambia.

The Trainers shall work with the Learners to collect evidence of competence, using the benchmarks provided by the unit standards. During the year, the Learners shall be required to undertake a series of practical assessment tasks. It is the sum of all these assessments tasks that deems a Learner to be competent (or not).

This approach to assessment is not a one-off event but one that gives learners many opportunities to demonstrate skill and allow for the capturing and recording of these demonstrations.

For the Learner to be deemed competent, they must demonstrate competency in every aspect of the practical tasks being undertaken. It must however be understood by the Trainer that Competency does not mean expert. It means that the candidate has attained sufficient skill and knowledge to perform the activity or service to a degree and quality that is acceptable to the industry and the customer in a time within which a competent person at the level could reasonably be expected to perform the task.

While this will be undertaken at institutional level, it is therefore envisaged that the Assessment principles of VALIDITY, RELIABILITY, FAIRENESS and FLEXIBILITY shall at all times be adhered to.



## Pre-Assessment

Assessment process explained to the employee (✓ if Yes).	<input type="checkbox"/>
Any appeal relating to the outcome of the assessment or the way in which the assessment was conducted shall be made through the company's <u>fair treatment policy</u> as explained to the employee (✓ if Yes).	<input type="checkbox"/>

<b>Employee/Trainee</b>  Employee/Trainee name: _____ (Print)  Employee/Trainee comments:		<b>Assessor</b>  Assessor name: _____ (Print)  Assessor comments:	
I fully understand the assessment and appeals process.		Theory assessment sighted and checked as satisfactory.	
Signature: _____  Date: _____		Signature: _____  Date: _____	

<b>Table of Contents</b>	<b>Page</b>
Technical Education, Vocational and Entrepreneurship Training Authority (TEVETA)	
(unit code) 671 Chemistry Techniques	1
<b>Record of Practical Assessment</b>	<b>1</b>
<b>Preface</b>	<b>2</b>
<b>Pre-assessment</b>	<b>3</b>
<b>Table of Contents</b>	<b>4</b>
<b>Prepare for practical assessment</b>	<b>5</b>
<b>Work health and safety</b>	<b>5</b>
<b>Customising the assessment</b>	<b>5</b>
<b>Carrying out the assessment</b>	<b>5</b>
<b>Completing the assessment</b>	<b>5</b>
<b>Assessor qualifications</b>	<b>5</b>
<b>Expiry status of assessment</b>	<b>5</b>
<b>Resources required</b>	<b>5</b>
<b>Range of variables</b>	<b>5</b>
<b>1. Safety in the Laboratory</b>	<b>6</b>
<b>2. Using Measuring Instruments</b>	<b>8</b>
<b>3. Determining the Water of Crystallisation</b>	<b>10</b>
<b>4. Using volumetric glassware</b>	<b>13</b>
<b>5. Weighing using an Analytical Balance</b>	<b>15</b>
<b>6. Preparation of a primary standard solution</b>	<b>18</b>
<b>7. Standardising a secondary standard</b>	<b>21</b>
<b>8. Determining the concentration of a sodium hydroxide solution</b>	<b>25</b>
<b>9. Sublimation of Ammonium Chloride</b>	<b>29</b>
<b>10. An experiment to show that a chemical reaction causes no change in mass</b>	<b>31</b>
<b>11. Separation of sand and sodium chloride</b>	<b>33</b>
<b>12. Final assessment summary</b>	<b>36</b>
<b>13. Assessment outcome</b>	<b>38</b>
<b>14. Validation of the assessment</b>	<b>39</b>



Prepare for the practical assessment

Add text here

Work Health and Safety

Add text here

Customising the assessment

Add text here

Carrying out the assessment

Add text here

Completing the assessment

Add text here

Assessor qualifications

Add text here

Expiry status of assessment

Add text here

Resources required

Add text here

Range of variables

Add text here

1.0 PERFORM AN EXPERIMENT ON SAFETY IN A LABORATORY	Satisfactory	Not Satisfactory
During observation of work activities, the candidate demonstrated that they can:		
a. Identify the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Dust mask</li> <li><input type="checkbox"/> Safety goggles</li> <li><input type="checkbox"/> Fume hood</li> <li><input type="checkbox"/> Latex gloves</li> <li><input type="checkbox"/> Gas mask (aspirators)</li> <li><input type="checkbox"/> Acid-resistant Laboratory coat</li> <li><input type="checkbox"/> Acid-resistant Safety shoes</li> <li><input type="checkbox"/> Fire extinguishers</li> <li><input type="checkbox"/> Fire blanket</li> <li><input type="checkbox"/> First-aid box</li> </ul> (No loose clothing, no long hair, no open shoes)	<input type="checkbox"/>	<input type="checkbox"/>
b. Set up apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Cleaning the identified pieces of apparatus</li> <li><input type="checkbox"/> Assembling the identified pieces of apparatus</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
c. Run the experiment correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Applying Dos and DONT's of the equipment and experiment</li> <li><input type="checkbox"/> Being alert of what is going on during the experiment</li> <li><input type="checkbox"/> Identifying the use of each equipment/apparatus</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
d. Interpret the results correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Appreciating the use and care of the equipment and apparatus</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
e. Identify the errors that may affect the experiment and how they can be avoided.	<input type="checkbox"/>	<input type="checkbox"/>



Assessor comments:

---

---

---

---

---

---

---

---

---

---

Signed:

Assessor:

---

Trainee:

---

2 PERFORMING AN EXPERIMENT ON USING MEASURING INSTRUMENTS	Satisfactory	Not Satisfactory
During observation of work activities, the candidate demonstrated that they can:		
a. Identify the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Analytical balance</li> <li><input type="checkbox"/> Top pan balance</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
b. Set up the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Placing the balances on the a strong and thick horizontal bench or concrete block</li> <li><input type="checkbox"/> Plugging in and switching on the balances</li> <li><input type="checkbox"/> Ensuring that the balance is on a horizontal surface by checking the bubble in the spirit level being in the centre of the circle.</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
c. Run the experiment correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Zeroing the apparatus when the sliding doors are closed</li> <li><input type="checkbox"/> Calibrating using the standard mass pieces</li> <li><input type="checkbox"/> Weighing out any material that is within the maximum load that the balance can withstand</li> <li><input type="checkbox"/> Weighing the mass either by taring or by difference</li> <li><input type="checkbox"/> Record the mass</li> <li><input type="checkbox"/> Repeating the process for four more materials and recording the masses</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
d. Interpret the results correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Checking the tolerance limits, deciding which is accurate or precise</li> <li><input type="checkbox"/> Checking for any bias in the measurement.</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
e. Identify the errors that may affect the experiment and how they can be avoided	<input type="checkbox"/>	<input type="checkbox"/>



Assessor comments:

---

---

---

---

---

---

---

---

Signed:

Assessor:

Trainee:

---

---

3 PERFORM AN EXPERIMENT TO DETERMINE THE WATER OF CRYSTALLISATION IN HYDRATED COPPER (II) SULPHATE	Satisfactory	Not Satisfactory
During observation of work activities, the candidate demonstrated that they can:		
a. Identify the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Clean dry porcelain crucible and lid</li> <li><input type="checkbox"/> Analytical balance</li> <li><input type="checkbox"/> Tripod stand and clay triangle or wire gauze</li> <li><input type="checkbox"/> Desiccator</li> <li><input type="checkbox"/> Clean crucible tongs</li> <li><input type="checkbox"/> Bunsen burner</li> <li><input type="checkbox"/> Gas cylinder</li> <li><input type="checkbox"/> Spatula</li> <li><input type="checkbox"/> Copper (II) sulphate pentahydrate</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
b. Set up the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Arranging the apparatus and materials in order</li> <li><input type="checkbox"/> Obtaining a clean crucible and lid.</li> <li><input type="checkbox"/> Checking the crucible for fractures. If none are found, support the crucible and lid on a clay triangle or wire gauze on a tripod stand</li> <li><input type="checkbox"/> Connecting the burner to the gas tap</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>



<p>c. Run the experiment correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Turning the gas on from the gas cylinder and lighting the Bunsen burner. Adjusting the air inlet vent until the flame is blue and the heat is intense</li> <li><input type="checkbox"/> Heating the crucible and lid for 5 minutes. Removing from the heat and placing in a desiccator to cool. Do not place on a laboratory bench to avoid contamination. Do not touch the crucible and lid as hot and cool crucibles look the same. Hence, handle them with crucible tongs. Handle the crucible and lid with the crucible tongs for the remainder of the experiment</li> <li><input type="checkbox"/> Determining the mass of the cool crucible and lid on an analytical balance and recording it, <math>m_1</math>.</li> <li><input type="checkbox"/> Weighing 3 g of hydrated copper (II) sulphate and transferring it to the crucible.</li> <li><input type="checkbox"/> Weighing and recording the combined mass of the crucible, lid and the hydrated copper (II) sulphate, <math>m_2</math>.</li> <li><input type="checkbox"/> Calculating the mass of the hydrated salt</li> <li><input type="checkbox"/> Returning the crucible (use crucible tongs only) with the sample to the clay triangle; set the lid just off the lip of the crucible to allow the evolved water molecules to escape on heating</li> <li><input type="checkbox"/> Heating the sample. Initially heating the sample slowly and then gradually intensifying the heat. Do not allow the crucible to become red hot. This may cause the anhydrous salt to decompose as well.</li> <li><input type="checkbox"/> Maintaining the high temperature on the sample for 10 minutes. Removing from the heat and placing in a desiccator to cool</li> <li><input type="checkbox"/> Determining the combined mass of the crucible, lid and anhydrous salt on the same balance used for earlier measurements, <math>m_3</math>.</li> <li><input type="checkbox"/> Determining if all the water has been removed by reheating the sample for 2 minutes, but do not intensify the flame</li> <li><input type="checkbox"/> Removing from the heat and placing in a desiccator to cool</li> <li><input type="checkbox"/> Measuring the combined mass again. If the second mass measurement of the anhydrous salt disagrees by greater than <math>\pm 0.010</math> g from the first, repeat the process of reheating the sample.</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
---	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------

<p>d. Interpret the results correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Calculating the percentage of the water of crystallization as follows:</li> <li><input type="checkbox"/> Mass of water of crystallization = <math>m_2 - m_3</math></li> <li><input type="checkbox"/> Mass of hydrated copper (II) sulphate = <math>m_2 - m_1</math></li> </ul> <p>% of water of crystallization = <math>\frac{m_2 - m_3}{m_2 - m_1} 100\%</math></p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
---	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------

Assessor comments:

---

---

---

---

---

---

---

---

---

---

Signed:

Assessor:

Trainee:

---



---



4PERFORMING AN EXPERIMENT ON HOW TO USE VOLUMETRIC GLASSWARE.	Satisfactory	Not Satisfactory
During observation of work activities, the candidate demonstrated that they can:		
a. Identifying the apparatus correctly.This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Pyrex</li> <li><input type="checkbox"/> Borosilicate</li> <li><input type="checkbox"/> Non-pyrex</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
b. Set up the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Cleaning with industrial detergent</li> <li><input type="checkbox"/> Washing with running tap water</li> <li><input type="checkbox"/> Rinsing with deionised water of distilled water</li> <li><input type="checkbox"/> Drying where necessary.</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
c. Run the experiment correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Filling the glassware with bigger volume requires transferring a liquid using a big beaker up to where the neck starts. Then topping up using a wash bottle to just slightly below the engraved graduation mark. Then finally using a jet to make up to the mark.</li> <li><input type="checkbox"/> Shaking well to homogenise the liquid</li> <li><input type="checkbox"/> Repeating the procedure for other types of volumetric glassware.</li> <li><input type="checkbox"/> Transferring the liquids into stoppered glass bottles</li> <li><input type="checkbox"/> Labelling includes:               <ul style="list-style-type: none"> <li>➤ Identity of the liquid</li> <li>➤ Identity of a person who prepared</li> <li>➤ Concentration of liquid where necessary</li> <li>➤ Date of preparation</li> <li>➤ Expiry date</li> </ul> </li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>





5 PERFORM AN EXPERIMENT ON WEIGHING USING AN ANALYTICAL BALANCE	Satisfactory			Not Satisfactory		
During observation of work activities, the candidate demonstrated that they can:						
a. Identify the apparatus/materials correctly. This may include:  <div><input type="checkbox"/> Analytical balance</div> <div><input type="checkbox"/> spatula,</div> <div><input type="checkbox"/> weighing boat or pan,</div> <div><input type="checkbox"/> beaker,</div> <div><input type="checkbox"/> stirring rod,</div> <div><input type="checkbox"/> suitable volumetric flask,</div> <div><input type="checkbox"/> distilled or deionised water,</div> <div><input type="checkbox"/> camel hair brush</div> <div><input type="checkbox"/> Primary standard eganalar anhydrous sodium carbonate,</div>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
b. Set up apparatus/ materials correctly. This may include:  <div><input type="checkbox"/> Calculating the mass of the primary standard material to be weighed based on relative atomic masses of the constituent atoms.</div> <div><input type="checkbox"/> Placing the balances on the a strong and thick horizontal bench or concrete block in a room where there are no air currents</div> <div><input type="checkbox"/> Plugging in the and switching on the balance</div> <div><input type="checkbox"/> Ensuring that the balance is on a horizontal surface by checking the bubble in the spirit level being in the centre of the circle.</div>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

<p>c. Run the experiment correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Placing the analytical balance on a horizontal surface by adjusting the screws on its base until the bubble on the spirit level is in the centre</li> <li><input type="checkbox"/> Weighing by taring includes <ul style="list-style-type: none"> <li>➤ Zeroing the balance,</li> <li>➤ weighing the empty receptacle,</li> <li>➤ Zeroing again</li> <li>➤ Putting the primary standard into the receptacle</li> <li>➤ Closing the sliding door</li> <li>➤ Waiting for the display to stabilize</li> <li>➤ Taking the reading</li> <li>➤ Recording the reading of the mass of the primary standard</li> </ul> </li> <li><input type="checkbox"/> Transferring the primary standard quantitatively into a clean dry beaker</li> <li><input type="checkbox"/> Dissolving it with minimum amount of distilled water</li> <li><input type="checkbox"/> Transferring the dissolved primary standard into a volumetric flask quantitatively</li> <li><input type="checkbox"/> Using distilled water from a clean beaker to add enough water until just below the neck of the volumetric flask</li> <li><input type="checkbox"/> Using a wash bottle to add more water to just below the engraved graduation mark.</li> <li><input type="checkbox"/> Using a jet to fill up to the graduation mark</li> <li><input type="checkbox"/> Shaking well to homogenise the solution</li> <li>Transferring the liquids into stoppered glass bottles</li> <li><input type="checkbox"/> Labelling includes: <ul style="list-style-type: none"> <li>➤ Identity of the liquid</li> <li>➤ Identity of a person who prepared</li> <li>➤ Concentration of liquid where necessary</li> <li>➤ Date of preparation</li> <li>➤ Expiry date</li> </ul> </li> </ul>						
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>



<p>d. Interpret the results correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Storing the prepared solutions in respective places.</li> <li><input type="checkbox"/> Checking the expiry date.</li> <li><input type="checkbox"/> Comparing the results from the solutions or liquids from time to time if necessary against secondary standard</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<p>e. Identify the errors that may affect the experiment and how they can be avoided</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Assessor comments:

---

---

---

---

---

---

---

---

---

---

Signed:      Assessor:                      Trainee:

\_\_\_\_\_

6 PERFORMING AN EXPERIMENT ON PREPARATION OF A STANDARD SOLUTION	Satisfactory	Not Satisfactory
During observation of work activities, the candidate demonstrated that they can:		
<p>a. Identify the apparatus /materials correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Beaker</li> <li><input type="checkbox"/> Stirring rod,</li> <li><input type="checkbox"/> Suitable</li> <li><input type="checkbox"/> Volumetric flask,</li> <li><input type="checkbox"/> Funnel.</li> <li><input type="checkbox"/> Primary standard materials eg anhydrous sodium carbonate,</li> <li><input type="checkbox"/> Distilled or deionised water</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
<p>b. Set up the apparatus correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Calculating the mass of the required standard from relative atomic masses.</li> <li><input type="checkbox"/> The weighed primary standard material (from the procedure on weighing) is transferred to another clean beaker where it is dissolved in the minimum amount of distilled or deionised water</li> <li><input type="checkbox"/> Selecting an appropriate volumetric flask and placing a clean ordinary funnel in its mouth with a piece of clean folded paper to serve as a wedge at its opening. This allows air to be easily expelled from the volumetric flask when the solution is poured into it.</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>



<p>c. Run the experiment correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Pouring the solution from the beaker into the funnel taking care not to spill any.</li> <li><input type="checkbox"/> Rinsing out the beaker with small portions of the water and adding to the volumetric flask each washing until everything has been washed out from the beaker.</li> <li><input type="checkbox"/> Rinsing, similarly, the funnel with several small amounts of water in a circular motion. Do not add such amounts so as to fill the volumetric flask.</li> <li><input type="checkbox"/> Removing the funnel from the volumetric flask and carefully and slowly adding enough distilled water to take the lowest point of the meniscus level with the engraved circular line on the neck of the flask.</li> <li><input type="checkbox"/> Using of a small jet of water is recommended when nearing the engraved line.</li> <li><input type="checkbox"/> Making an allowance for water from above the line especially on the inside of the neck to roll down as ignoring this may lead to overshooting.</li> <li><input type="checkbox"/> Stoppering the volumetric flask and shaking vigorously in order to prepare a homogeneous solution.</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<p>d. Interpret the results correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Titrating against secondary standards</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<p>e. Identify the errors that may affect the experiment and how they may be avoided</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Assessor comments:

---

---

---

---

---

---

---

---

---

---

Signed:

Assessor:

Trainee:

---

---



7. PERFORMING AN EXPERIMENT TO STANDARDISE A SECONDARY STANDARD AGAINST A PRIMARY SRANDARD.	Satisfactory			Not Satisfactory		
During observation of work activities, the candidate demonstrated that they can:						
a. Identify the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Pipette</li> <li><input type="checkbox"/> Burette</li> <li><input type="checkbox"/> Conical flasks</li> <li><input type="checkbox"/> Volumetric flask</li> <li><input type="checkbox"/> Measuring cylinder</li> <li><input type="checkbox"/> Pipette filler</li> <li><input type="checkbox"/> Beakers including 1 litre beaker for diluting</li> <li><input type="checkbox"/> Clamp and stand</li> <li><input type="checkbox"/> White tile</li> <li><input type="checkbox"/> Concentrated sulphuric acid (from which to prepare the secondary standard)</li> <li><input type="checkbox"/> Sodium carbonate standard solution – primary standard</li> <li><input type="checkbox"/> Phenolphthalein indicator (other suitable indicators like methyl orange or methyl red may be used)</li> <li><input type="checkbox"/> Distilled water</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

<p>b. Set up the experiment correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Cleaning apparatus using industrial detergent</li> <li><input type="checkbox"/> Rinsing apparatus using running tap water</li> <li><input type="checkbox"/> Rinsing apparatus using distilled water</li> <li><input type="checkbox"/> Drying apparatus where necessary</li> <li><input type="checkbox"/> Assembling apparatus</li> <li><input type="checkbox"/> Most concentrated acid bottles have the following data on their labels;</li> <li><input type="checkbox"/> Percentage purity, specific gravity SG, and molecular weight of that particular acid.</li> <li><input type="checkbox"/> Thus, the following formula is used to calculate the concentration of the concentrated acid.</li> <li><input type="checkbox"/> Molarity = %purity XSG X1000/molecular weight</li> <li><input type="checkbox"/> Using the dilution <math>M_1V_1 = M_2V_2</math> where the left side refers to the concentrated solution while the right side pertains to the dilute solution. M represents the molarities while V represents the volumes.</li> <li><input type="checkbox"/> <b>Note to the lecturer: dilution of sulphuric acid is quite dangerous.</b></li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
--	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------	--------------------------



<p>c. Run the experiment correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Using a pipette with a pipette filler to draw three 25 ml aliquots from the dilute solution and placing them in clean, dry separate conical flasks. Adding 2 to 3 drops of phenolphthalein indicator to each of the conical flasks.</li> <li><input type="checkbox"/> Filling a clean burette with the sodium carbonate solution from tip to top.</li> <li><input type="checkbox"/> Adjusting the volume to 0.00 ml on the burette. Placing a white tile or white piece of paper below the burette on the base of the clamp.</li> <li><input type="checkbox"/> Placing a conical flask below the tip of the burette and continuing titrating until a persistent faint pink colour appears in the conical flask. The stopcock should be closed immediately at this point and the volume of the sodium carbonate solution dispensed from the burette read off immediately. The eyes should be at 90° to the burette at the point where the meniscus is. Repeating the procedure with the remaining two conical flasks. The volume readings should be entered in a table as shown below and they should agree to within 0.05 ml of each other.</li> <li><input type="checkbox"/> Running, slowly, the sodium carbonate solution into the conical flask while swirling all the time. The right hand should handle the conical flask while the left hand should operate the stopcock</li> </ul>																									
<p>d. Interpret the results correctly. This may include:</p> <table border="1" data-bbox="199 1391 962 1731"> <thead> <tr> <th></th><th>Initial volume, ml</th><th>Final volume, ml</th><th>Actual volume, ml</th></tr> </thead> <tbody> <tr> <td>Trial 1</td><td></td><td></td><td></td></tr> <tr> <td>Trial 2</td><td></td><td></td><td></td></tr> <tr> <td>Trial 3</td><td></td><td></td><td></td></tr> </tbody> </table> <p><input type="checkbox"/> Use the following chemical equation for the calculations:</p> $\text{Na}_2\text{CO}_{3(aq)} + \text{H}_2\text{SO}_{4(aq)} \rightarrow \text{Na}_2\text{SO}_{4(aq)} + \text{H}_2\text{O}_{(l)} + \text{CO}_{2(g)}$					Initial volume, ml	Final volume, ml	Actual volume, ml	Trial 1				Trial 2				Trial 3									
	Initial volume, ml	Final volume, ml	Actual volume, ml																						
Trial 1																									
Trial 2																									
Trial 3																									

Assessor comments:

---

---

---

---

---

---

---

---

---

---

Signed:

Assessor:

Trainee:

---

---

25

<p>b. Set up apparatus correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Cleaning apparatus using industrial detergent</li> <li><input type="checkbox"/> Rinsing apparatus using running tap water</li> <li><input type="checkbox"/> Rinsing apparatus using distilled water</li> <li><input type="checkbox"/> Drying apparatus where necessary</li> <li><input type="checkbox"/> Assembling apparatus</li> <li><input type="checkbox"/> Most concentrated acid bottles have the following data on their label</li> <li><input type="checkbox"/> Percentage purity, specific gravity SG, and molecular weight of that particular acid.</li> <li><input type="checkbox"/> Thus, the following formula is used to calculate the concentration of the concentrated acid.</li> <li><input type="checkbox"/> Molarity = %purity XG X1000/molecular weight <ul style="list-style-type: none"> <li><input type="checkbox"/> Using the dilution <math>M_1V_1 = M_2V_2</math> where the left side refer to the concentrated solution while the right side pertain to the dilute solution. M represents the molarities while V represents the volumes.</li> <li><input type="checkbox"/> Removing <math>V_1</math>, that is, the calculated volume of the concentrated acid that should be taken and diluted. Adding the acid to the glass funnel in small quantities with constant stirring with a glass rod before adding the next instalment. Continuing stirring for some time and leaving to cool.</li> <li><input type="checkbox"/> Transferring to a 1 litre volumetric flask and making up to volume with distilled water.</li> <li><input type="checkbox"/> Shaking the flask well to homogenise the solution.</li> <li><input type="checkbox"/> The concentration of the final solution should be 0.50 M. This is just a guide as other concentration may be prepared depending on the requirement</li> <li><input type="checkbox"/> Preparing the sodium hydroxide solution:</li> <li><input type="checkbox"/> Weighing about 20 g of sodium hydroxide on a watch glass or small clean dry beaker quickly as sodium hydroxide absorbs moisture from the atmosphere.</li> <li><input type="checkbox"/> Transferring the pellets to another clean beaker and dissolving in the minimum amount of distilled water.</li> <li><input type="checkbox"/> Transferring the sodium hydroxide solution to a clean 1 litre volumetric flask and adding enough distilled water to make up to the mark.</li> <li><input type="checkbox"/> Shake the volumetric flask well to homogenise.</li> </ul> </li> </ul>						
--	--	--	--	--	--	--



<p>c. Run the experiment correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Using a pipette with a pipette filler to draw three 25 ml aliquots from the dilute solution of sodium hydroxide and placing them in clean, dry separate conical flasks. Adding 2 to 3 drops of phenolphthalein indicator to each of the conical flasks. The colour of the solution will change to pink.</li> <li><input type="checkbox"/> Filling a clean burette with the hydrochloric acid standard solution from tip to top.</li> <li><input type="checkbox"/> Adjusting the volume to 0.00 ml on the burette. Placing a white tile or white piece of paper below the burette on the base of the clamp.</li> <li><input type="checkbox"/> Placing a conical flask below the tip of the burette and slowly running the hydrochloric acid solution into the conical flask while swirling all the time. The right hand should handle the conical flask while the left hand should operate the stopcock.</li> <li><input type="checkbox"/> Continuing titrating until there is a distinct colour change in the conical flask, i.e the solution becomes colourless. The stop cork should be closed immediately at this point and the volume of the hydrochloric acid solution dispensed from the burette read-off immediately. The eyes should be at 90° to the burette at the bottom of the meniscus.</li> <li><input type="checkbox"/> Repeating the procedure with the remaining two conical flasks. The volume reading should be entered in a table as shown below and they should agree to within 0.05 ml of each other</li> <li><input type="checkbox"/> Note that due to advances there exists autotitrators and other means of indicating the end other than indicators.</li> </ul>																						
<table border="1" data-bbox="284 1512 963 1879"> <thead> <tr> <th></th><th>Initial volume, ml</th><th>Final volume, ml</th><th>Actual volume, ml</th></tr> </thead> <tbody> <tr> <td>Trial 1</td><td></td><td></td><td></td></tr> <tr> <td>Trial 2</td><td></td><td></td><td></td></tr> <tr> <td>Trial 3</td><td></td><td></td><td></td></tr> </tbody> </table> <p>Copy and complete the table</p>		Initial volume, ml	Final volume, ml	Actual volume, ml	Trial 1				Trial 2				Trial 3				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Initial volume, ml	Final volume, ml	Actual volume, ml																			
Trial 1																						
Trial 2																						
Trial 3																						

d. Interpret the results correctly. This may include:  Use the equation:  $\text{NaOH}_{(\text{aq})} + \text{HCl}_{(\text{aq})} \rightarrow \text{NaCl}_{(\text{aq})} + \text{H}_2\text{O}_{(\text{l})}$  <input type="checkbox"/> Examining the individual volumes for precision. <input type="checkbox"/> Calculating individual results for each volume <input type="checkbox"/> Comparing the individual results for precision	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
e. Identify the error that may affect the experiment and how they may be avoided	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Assessor comments:

---

---

---

---

---

---

---

---

---

---

Signed:

Assessor:

Trainee:

---



---



9PERFORMING AN EXPERIMENT ON THE SUBLIMATION OF AMMONIUM CHLORIDE	Satisfactory	Not Satisfactory
During observation of work activities, the candidate demonstrated that they can:		
a. Identify the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Pyrex test-tube</li> <li><input type="checkbox"/> Analytical balance</li> <li><input type="checkbox"/> Small clean beaker</li> <li><input type="checkbox"/> Bunsen burner or spirit lamp</li> <li><input type="checkbox"/> Clamp and stand</li> <li><input type="checkbox"/> Gas cylinder:</li> <li><input type="checkbox"/> Anhydrous ammonium chloride powder</li> <li><input type="checkbox"/> Box of matches</li> <li><input type="checkbox"/> Sand</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
b. Set up the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Weighing about 0.5 g of ammonium chloride in a small clean and dry beaker</li> <li><input type="checkbox"/> Transferring quantitative the ammonium chloride to a clean dry test-tube</li> <li><input type="checkbox"/> Clamping the test-tube above the Bunsen burner</li> <li><input type="checkbox"/> Turning the gas on from the gas cylinder</li> <li><input type="checkbox"/> Bringing the lit match stick close to the top of the Bunsen burner to light the gas</li> <li><input type="checkbox"/> Adjusting the air inlet vent at the base of the burner until you have a bluish flame</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
c. Run the experiment correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Heating gently the bottom of the test-tube. The ammonium chloride will disappear completely from the bottom of the test-tube after a while. Some escapes as vapour but condenses on the upper wall of the test-tube.</li> <li><input type="checkbox"/> Repeating the experiment with the ammonium chloride mixed with 1 g of sand and noting your observations.</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
d. Interpret the results correctly. This may include:	<input type="checkbox"/>	<input type="checkbox"/>
e. Identify the errors that may have affected the experiment and how they can be avoided.	<input type="checkbox"/>	<input type="checkbox"/>

Assessor comments:

---

---

---

---

---

---

---

---

---

---

Signed:	Assessor:	Trainee:
	<hr/>	<hr/>



10 PERFORMING AN EXPERIMENT TO SHOW THAT A CHEMICAL REACTION CAUSES NO CHANGE IN WEIGHT	Satisfactory	Not Satisfactory
During observation of work activities, the candidate demonstrated that they can:		
a. Identify the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Analytical balance</li> <li><input type="checkbox"/> Small clean test-tube</li> <li><input type="checkbox"/> Clean 600 ml beaker</li> <li><input type="checkbox"/> Flexible thin wire</li> <li><input type="checkbox"/> Small clean beaker</li> <li><input type="checkbox"/> Barium chloride or barium nitrate</li> <li><input type="checkbox"/> Sodium sulphate</li> <li><input type="checkbox"/> Distilled or deionised water</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
b. Set up the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Cleaning some apparatus industrial detergent</li> <li><input type="checkbox"/> Rinsing some apparatus with running tap water</li> <li><input type="checkbox"/> Rinsing some apparatus with distilled water</li> <li><input type="checkbox"/> Assembling some apparatus</li> <li><input type="checkbox"/> Weighing</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>

<p>c. Running the experiment. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Weighing the 600 ml beaker</li> <li><input type="checkbox"/> Adding about 3 g of barium chloride to the beaker and record the weight</li> <li><input type="checkbox"/> Weighing the small test-tube</li> <li><input type="checkbox"/> Weighing about 3 g of sodium sulphate and transferring it to the small test-tube</li> <li><input type="checkbox"/> Weighing the test-tube again and record the new weight</li> <li><input type="checkbox"/> Dissolving the barium chloride in the minimum quantity of water. Weighing the beaker again and recording the weight</li> <li><input type="checkbox"/> Dissolving the sodium sulphate in the minimum amount of water</li> <li><input type="checkbox"/> Tying the flexible thin wire around the upper part of the test-tube</li> <li><input type="checkbox"/> Placing the test-tube in the beaker in an inclined position taking care not to bring the two solutions into contact. Hooking one end of the wire on the rim of the 600 ml beaker</li> <li><input type="checkbox"/> Weighing carefully this assembly again and recording the mass, <math>m_1</math>.</li> <li><input type="checkbox"/> Removing the assembly from the balance</li> <li><input type="checkbox"/> Unhooking the wire and letting all the contents of the test-tube to come into contact with the barium chloride solution in the beaker. A white precipitate will form. Wait for 10 minutes and reweigh the beaker and everything in it again, <math>m_2</math>.</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<p>d. Interpret data correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Comparing <math>m_1</math> and <math>m_2</math></li> <li><input type="checkbox"/> <math>m_1</math> and <math>m_2</math> should agree</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<p>e. Identify errors that may affected the experiment and how they can be avoided</p>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Assessor comments:

---



---



---



---

Signed:      Assessor:      Trainee:

\_\_\_\_\_



11 PERFORMING AN EXPERIMENT ON THE SEPARATION OF SAND AND SODIUM CHLORIDE.	Satisfactory	Not Satisfactory
During observation of work activities, the candidate demonstrated that they can:		
a. Identify the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Analytical balance</li> <li><input type="checkbox"/> Conical flask</li> <li><input type="checkbox"/> Evaporating dish</li> <li><input type="checkbox"/> Ordinary funnel</li> <li><input type="checkbox"/> Beakers</li> <li><input type="checkbox"/> Filter paper</li> <li><input type="checkbox"/> Clean evaporating dish (optional)</li> <li><input type="checkbox"/> Spirit lamp (optional)</li> <li><input type="checkbox"/> Tripod stand and wire gauze (optional)</li> <li><input type="checkbox"/> Sodium chloride powder</li> <li><input type="checkbox"/> Finely divided sand</li> <li><input type="checkbox"/> Distilled or deionised water</li> <li><input type="checkbox"/> oven may be used for speedy drying.</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
b. Set up the apparatus correctly. This may include: <ul style="list-style-type: none"> <li><input type="checkbox"/> Washing some apparatus with industrial detergent</li> <li><input type="checkbox"/> Rinsing some apparatus with running tap water</li> <li><input type="checkbox"/> Rinsing some apparatus with distilled water</li> <li><input type="checkbox"/> Drying some apparatus</li> <li><input type="checkbox"/></li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>

<p>c. Run the experiment correctly. This may include:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Weighing about 3 g of sodium chloride in a beaker</li> <li><input type="checkbox"/> Weighing 3 g of sand and it to beaker with the sodium chloride solution.</li> <li><input type="checkbox"/> Adding more water to the beaker until half way up.</li> <li><input type="checkbox"/> Putting a clean funnel on the conical flask with a small clean folded piece of paper to serve as a wedge for easy expulsion of air from the conical flask.</li> <li><input type="checkbox"/> Folding a clean filter paper and place it on the funnel with the cone resting at the base. The filter paper should be of such quality that it will retain the particles of sand. Agitating the sodium chloride/sand mixture.</li> <li><input type="checkbox"/> Pouring carefully out the mixture onto the filter paper. The sand will tend to remain in the beaker so agitate every time before pouring onto to the filter paper. Whatever remains in the beaker should be washed out with small quantities of water until everything has been transferred from the beaker.</li> <li><input type="checkbox"/> If solid sodium chloride is required, transfer the filtrate in the conical flask to an evaporating dish and heat on wire gauze to drive out most of the water. Do not heat to dryness. Remove the evaporating dish from the heat and leave it in a cool place where it won't be disturbed for 4 hours. After cooling, crystals will form. Decant any excess water and transfer the crystals to another clean filter paper and leave to dry. If dry sand is required, transfer the residue to a watch glass or filter paper and leave to dry. An oven may be used for speedy drying</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
d. Interpret the data correctly. This may include:	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
e. Identify errors that may have affected the experiment and how they can be avoided.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>



Assessor comments:

---

---

---

---

---

---

---

---

---

---

Signed:      Assessor:                      Trainee:

## Final Assessment Summary

### Practical assessment summary

Note: refer to mapping document if required

		Satisfactory	Not Satisfactory
1.	Safety in the Laboratory	<input type="checkbox"/>	<input type="checkbox"/>
2.	Measurements	<input type="checkbox"/>	<input type="checkbox"/>
3.	Determine The Water of Crystallisation Of Hydrated Copper (II) Sulphate	<input type="checkbox"/>	<input type="checkbox"/>
4.	Using Volumetric Glassware	<input type="checkbox"/>	<input type="checkbox"/>
5.	Weighing On An Analytical Balance	<input type="checkbox"/>	<input type="checkbox"/>
6.	Preparation Of A Primary Standard Solution	<input type="checkbox"/>	<input type="checkbox"/>
7.	Standardising A Secondary Standard	<input type="checkbox"/>	<input type="checkbox"/>
8.	Determining The Concentration Of Sodium Hydroxide	<input type="checkbox"/>	<input type="checkbox"/>
9.	Sublimation Of Ammonium Chloride	<input type="checkbox"/>	<input type="checkbox"/>
10.	A Chemical Reaction Causes No Change In Mass	<input type="checkbox"/>	<input type="checkbox"/>
11.	Separation Of Sand And Sodium Chloride	<input type="checkbox"/>	<input type="checkbox"/>

[illegible]

Signed:      Assessor:                      Trainee:

\_\_\_\_\_

## Assessment Outcome

Satisfactory ☐

Not Satisfactory ☐

<b>Employee/Trainee</b>	<b>Assessor</b>
Employee/Trainee name: _____ (Print)	Assessor name: _____ (Print)
<b>Employee/Trainee comments:</b>	<b>Assessor comments:</b>
Signature: _____  Date: _____	Signature: _____  Date: _____



## VALIDATION OF THE ASSESSMENT

NAME:..... DATE:.....

POSITION: **PRINCIPAL/HEAD OF INSTITUTION** SIGNATURE:.....

NAME INSTITUTION:.....

STAMP:

